

# Acetic Acid 5% & 8%

MenidiMedica  
Biotechnology Applications

**REF** 70901-100: 100 mL.  
70901-200: 200 mL.  
70901-500: 500 mL.  
70901-1000: 1 lt.



24 months, storage at RT



## Περιεχόμενο

Οξικό οξύ, Συντηρητικά

## Περιγραφή

Οξικό οξύ

## Τρόπος διάθεσης

Γυάλινες φιάλες των 100, 200, 500, 1000 mL.

## Δράση & χρήση

Το ACETIC ACID είναι ένα μέσο διαφοροποίησης των φυσιολογικών από τους επίφοβους ιστούς. Οι φυσιολογικοί ιστοί παραμένουν ως έχουν, ενώ ο ιστός που είναι ύποπτος για δυσπλασίες αναπτύσσει περιοχές λευκές. Μετά τη χρήση, πλύνετε με νερό ή φυσιολογικό ορό.

## Προειδοποιήσεις

**Μόνο για Εξωτερική Χρήση.** Σε περίπτωση ακούσιας οφθαλμικής επαφής, τα μάτια θα πρέπει να πλυθούν αμέσως με άφθονο νερό ή φυσιολογικό ορό, σηκώνοντας κατά διαστήματα τα άνω και κάτω βλέφαρα μέχρις ότου δεν υπάρχει καμία εμφάνιση υπολειμμάτων του διαλύματος (περίπου 15-20 λεπτά).

## Αποθήκευση

Διατηρείτε το καλά κλεισμένο και προστατευμένο από το φως. ΜΗΝ το χρησιμοποιείτε εάν η σφραγίδα ασφαλείας έχει σχιστεί. Αποθηκεύστε το σε θερμοκρασία δωματίου 15°-30°C.

## Απόρριψη

Οι ανοιγμένες φιάλες με αχρησιμοποίητο μέρος του προϊόντος που περιέχουν υπολείμματα του προϊόντος θα πρέπει να τοποθετούνται σε ένα κατάλληλο, ξηρό περιέκτη, ώστε ν'απορριφθούν αργότερα σύμφωνα με τις τοπικές πρακτικές διαχείρισης επικίνδυνων αποβλήτων.



## Content

Acetic Acid, Preservative

## Description

Acetic Acid

## How Supplied

Glass bottles of 100, 200, 500, 1000 mL.

## Action & Uses

Acetic acid is an agent used for differentiating normal from suspicious tissues. Normal tissue remain as they are, while tissue suspicious for malformations develops areas of white. Wash after use with water or normal saline.

## Warnings

**For External Use Only.** Should inadvertent ocular administration occur, the eye(s) should be washed immediately with large amounts of water or normal saline, occasionally lifting the upper and lower lids until no evidence of solution remains (approximately 15-20 minutes).

## Storage

Keep tightly closed and protect from light. DO NOT use if seal is broken. Store at room temperature 15°- 30°C.

## Disposal

Opened containers with unused portions of product should be placed in a suitable, dry container for later disposal according to local hazardous waste practices.



MenidiMedica Biotech  
Menidi Aetolias-Akarnanias  
30016-Greece  
+30 2681088000, +30 6937115868  
menidimedica@gmail.com  
www.menidimedica.g

